

## EMDIRENAL CL

***General purpose lightfast polymer with high dispersing and stabilizing power to improve grain tightness, softness and fullness.***

### **Descriptive properties**

<b>Composition</b>	: Addition polymerisation product
<b>Appearance</b>	: Brownish clear liquid
<b>pH (1:10)</b>	: 5,0 - 6,0
<b>Active content</b>	: App. 40 %
<b>Solubility in water</b>	: Easily soluble in water at 35°C
<b>Storage</b>	: Minimum 1 year - dry storage-

Emdirenal CL has a wide range of applications but its most important characteristic is its stabilizing power: Emdirenal CL improves softness and tightness, accelerates the distribution and penetration of anionic chemicals through the leather cross-section, improves dry-drumming resistance of the grain and promotes very good dyeing properties.

### **Application recipes**

In the beginning of retanning, based on wetblue shaved weight:

0,5 - 2,0 % Emdirenal CL

1,5 - 6,0 % Acrylic copolymer

### **Remarks**

Please refer to its MSDS.

### **Packaging**

120 kg in polypropylene drums.