

EMDIRENAL FIXT

Strong fixing auxiliary to yield full shades with improved wet and perspiration fastness with anionic dyes.

Descriptive properties

| | |
|------------------------------------|-----------------------------------|
| Appearance | : Fine white powder |
| pH (1:10) | : App. 4,2 |
| Al₂O₃ | : % 47,4 ± 1,0 |
| Ionicity | : Strongly cationic |
| Solubility in water | : Easily soluble in water at 35°C |
| Storage | : Minimum 1 year - dry storage- |

Emdirenal FIXT has a wide range of applications but its most important characteristic is its distinctive power for fixing anionic products: it increases the affinity of anionic products to the leather.

Emdirenal FIXT improves dye fastness and deepens the colour shade.

With the use of Emdirenal FIXT a darker and more even colour shade out of vegetable or synthetic retanned wetblue is obtained. The touch is improved.

Emdirenal FIXT can be used in fixing of waterproof fatliquors for white leathers.

Application recipes

Emdirenal FIXT should be diluted 1:6 - 1:8 with water prior to its addition to the float.

To improve general fastness:

In dyeing of aniline leathers: 0,2 - 1,0 % based on shaved weight

In dyeing of suede leathers: 1,0 - 3,0 % based on crust weight

In fixation of waterproof leathers: 1,5 - 2,0 % based on shaved weight

Remarks

Always use in fresh bath.

Packaging

20 kg in kraft bags.